





Recare 
China Tianjin Recare Co., Ltd.




Company Introduction

China Tianjin Recare Co.Ltd. is a pioneer in the condom manufacturing industry and distribution in China. Condoms supplied by Recare are highly regarded for the its durability and reliability in world wide government tenders, commercial markets, non commercial marketo Recare know only quality condoms provide reliable protection, in order to fulfil the highest quality, we equipped our factories with a team of experienced and qualified process engineers and technologists, to continuously improve on machinery,manufacturing processes. We have been assessed and certified to the requirement of ISO 9001,ISO13485,Directive 93/42/EEC.

All prodction processes are monitored permanently, data are regularly collected and evaluated. A certain number of condoms from each batch are selected to undergo further tests. Condoms from each batch are “archived”. With the help of the batch nubmer, it is possible to determine the time of manufacture, the machines responsible, and the materials used for each individual condom.

So we promised HIGH QUALITY, PERFECT SERVICE, PROMPT DELIVERYTIME, COMPETITIVE PRICE. We are your best supplier.

CONTRACEPTION

Types	FOIL POUCH SIZE (mm)	LENGTH (mm)	WIDTH (mm)	THICKNESS (mm)	COLOR	LUBRICANT (mg)
 classic	55X55 or 56x58	180 185 190 200	46±2 49±2 53±2 55±2 56±2	0.04±0.015 0.05±0.015 0.06±0.015 0.07±0.015	natural	≥250 ≥450 ≥650 ≥850
 ribbed	55X55 or 56x58	180 185	52±2 55±2	0.06±0.015	natural	Silicone Oil≥ 350
 dotted	55X55 or 56x58	180 185	49±2 52±2 55±2	0.06±0.015	natural	Silicone Oil≥ 350
 Contoured	55X55 or 56x58	180 185	49±2 52±2	0.06±0.015	natural	Silicone Oil≥ 350
 color	55X55 or 56x58	180 185	52±2	0.06±0.015	natural red yellow blue	Silicone Oil≥ 350
 delay	55X55 or 56x58	180 185	52±2	0.06±0.015	natural	Silicone Oil≥ 350
 ultra thin	55X55 or 56x58	180 185	52±2	0.05 0.045 0.035	natural	Silicone Oil≥ 350
 ultra big	55X55 or 56x58	185 190 200	56±2	0.06±0.015	natural	Silicone Oil≥ 350
 Graphene	55X55 or 56x58	180 - 185	52±2	0.03	natural	Silicone Oil≥ 350

Recare 001 Polyurethane Condom

- Packing: 3 pcs / box/12 pcs / box/24 pcs / box
- Description: ultra thin polyurethane condom
- Width: 55+/-2mm
- Length: 175-180mm
- Thickness: 0.01mm、0.02mm、0.03mm
- Water-soluble polyurethane
- Certificate: CE ISO FREE SALE



Recare Graphene condom

- Packing: 3 pcs / box
- Description: Latex and Graphene
- Width: 52+/-2mm
- Length: 185+/-2mm
- Silicone oil
- thickness: 0.035-0.04



Recare Classic Condom

- Packing: 3pieces/box, 12pieces/box
- Description: regular shape condom
- Width: 52+-2mm, 54+-2mm
- Length: 170mm, 180mm, 185mm
- Rectangular or Square foil pouch as you desire
- Spermicidal type available if needed
- Standard : ISO4074 2002.



Recare Ribbed Condom

- Packing: 3pieces/box, 12pieces/box
- Description: Textured for extra sensation
- Width: 52+-2mm, 54+-2mm
- Length: 170mm, 180mm, 185mm
- Rectangular or Square foil pouch as you desire
- Spermicidal type available if needed
- Standard : ISO4074 2002.



Recare Dotted Condom

- Packing: 3pieces/box, 12pieces/box
- Description: Dotted for extra sensation
- Width: 52+-2mm, 54+-2mm
- Length: 170mm, 180mm, 185mm
- Rectangular or Square foil pouch as you desire
- Spermicidal type available if needed
- Standard : ISO4074 2002.



Recare Fruit & Colour Condom

- ◎ Packing: 3pieces/box,12pieces/box
- ◎ Description: Fruit scented and coloured latex
- ◎ Width:52+-2mm,54+-2mm
- ◎ Length:170mm,180mm,185mm
- ◎ Rectangular or Square foil pouch as you desire
- ◎ Spermicidal type available if needed
- ◎ Standard : ISO4074 2002.



Recare Contoured Condom(Ribbed and Dotted Condom)

- ◎ Packing:3pieces/box,12pieces/box
- ◎ Description:dotted and ribbed for extra sensation
- ◎ Width:52+-2mm,54+-2mm
- ◎ Length:170mm,180mm,185mm
- ◎ Rectangular or Square foil pouch as you desire
- ◎ Spermicidal type available if needed
- ◎ Standard : ISO4074 2002.



Recare Delay Condom

- ◎ Packing: 3pieces/box,12pieces/box
- ◎ Description:Contains an active ingredient called benzocaine, which make mildly and Temporarily desensitizes the nerve endings in the penis,prolonging Love making.
- ◎ Width:52+-2mm,54+-2mm
- ◎ Length:170mm,180mm,185mm
- ◎ Rectangular or Square foil pouch as you desire
- ◎ Standard : ISO4074 2002.



Recare Ultra Thin Condom

- ◎ Packing: 3pieces/box,12pieces/box
- ◎ Description: Ultra thin condoms are ultra thin for a more natural feeling. Ideal for the experience user, seeking extra sensation, they Are thinner than all other condoms, almost like using nothing at all.
- ◎ Width:52+-2mm,54+-2mm
- ◎ Length:170mm,180mm,185mm
- ◎ Rectangular or Square foil pouch as you desire
- ◎ Standard : ISO4074 2002.



Recare Ultra Big Condom

- ◎ Packing: 3pieces/box,6pieces/box,12pieces/box
- ◎ Description: Latex condoms are intended for men who feel that current Regular and large size condoms are too small, Other men may Experience slippage with this extra large size condom
- ◎ Width: 56+-2mm
- ◎ Length: 190mm-195mm
- ◎ Square foil pouch
- ◎ Standard : ISO4074 2002.



Recare Ultra Thick Condom

- ◎ Packing: 3pieces/box,6pieces/box,12pieces/box
- ◎ Description: Ultra-Thick Latex Condom is designed for people who like to practice safe sex. The extra thickness on the supreme quality latex can help to give extra barrier from many sexually transmitted disease, yet maintaining sensitivity. It is developed for sustained enjoyment in mind.
- ◎ Width: 54+-2mm, 56+-2mm
- ◎ Length: 185mm,190mm ,190mm-195mm
- ◎ thickness: 0.07mm
- ◎ Square foil pouch
- ◎ Standard : ISO4074 2002.



Recare Freezing feeling Condom

- ◎ Packing: 3pieces/wallet,3pieces/box,6pieces/box,12pieces/box
- ◎ Description: These condoms are lubricated in a freezing element,create a unprecedented comfortable feeling. That provide you the extra pleasure.
- ◎ Width: 54+-2mm, 56+-2mm
- ◎ Length: 185mm,190mm
- ◎ Rectangular or Square foil pouch as you desire
- ◎ Standard : ISO4074 2002.



Recare Warming feeling Condom

- ◎ Packing: 3pieces/wallet,3pieces/box,6pieces/box,12pieces/box
- ◎ Description: These condoms are lubricated in a warming element,create a unprecedented comfortable feeling. That provide you the extra pleasure.
- ◎ Width: 54+-2mm, 56+-2mm
- ◎ Length: 180mm,185mm,190mm
- ◎ Rectangular or Square foil pouch as you desire
- ◎ Standard : ISO4074 2002.



Recare Ultra Lubricant Condom

- ◎ Packing: 3pieces/box,6pieces/box,12pieces/box
- ◎ Description: More lubricated to overcome friction and created a great natural feeling
- ◎ Width: 52+-2mm,54+-2mm
- ◎ Length: 180mm,185mm,190mm
- ◎ Square foil pouch
- ◎ Standard : ISO4074 2002.



Recare spike condom

- ◎ Packing: 1pieces/box
- ◎ Description: Spike condom
- ◎ Width: 52+-2mm,54+-2mm
- ◎ Length: 172-192mm
- ◎ Pink,blue,green,yellow,red,black or as required
- ◎ Certification:CE,ISO,FSC,SABS



Recare vegan condom

- ◎ Processing plant custom made free casein green vegan
- ◎ Description: Sex Products
- ◎ Width: 52+-2mm,54+-2mm
- ◎ Length: 180-192mm
- ◎ Natural / Pink / Blue / Green / Yellow etc
- ◎ Certification:CE,ISO,FSC,SABS



Recare Finger condom

- ◎ Professional factory oem service higher quality finger condom
- ◎ Description: Natural latex
- ◎ Width: 25-29mm
- ◎ Length: 75mm-85mm
- ◎ Natural / Pink / Blue / Green / Yellow etc
- ◎ CE, ISO FSC, SABS, SGS, Directive 93/42/EEC



INTRAUTERINE DEVICE (IUD) Kit



IUD Kit
Product Information 022365 IUD Kit
(25/Box)
022365-M Metal Disposable IUD Kit
(25/Box)

COPPERT

The TCu 380A intrauterine device (IUD) is one of the most effective, long-acting reversible contraceptives available. It has more copper surface (380 square millimetres) than any other marketed IUD and is the most advanced of the family of T-shaped copper devices developed by the population council.

The TCu 380A IUD is inserted high in uterus by a physician, providing enhanced contraceptive for up to ten years. The annual pregnancy rate of users of the TCu 380A is well below one per 100 users.

The TCu 380A IUD consists of flexible polyethylene T-shaped device with a copper collar of approximately 68.7 mg of copper on each of its transverse arms and 176 mg of copper wire coiled around its vertical shaft. A monofilament polyethylene thread is tied through the bulb at the tip of the T, resulting in two threads to aid in detection and removal of the device.





FOR ORDERS
Email: supplies@kumis.co.za
www.kumis.co.za