

DENTAL





HS-M01 Dental Bib

Material: 2x20g Paper Laminated 20g Pe Film
Size: 38x45cm, 33x66cm or 40x45cm
Color: white, green, pink, blue, green, purple
Packing: 100pcs/bag, 500pcs/ctn
Carton Size: 34x29x25cm



HS-M02 Saliva Ejector


Material: Medical Grade PVC
Length: 15cm
Packing: 100pcs/bag, 10bags/ctn
Carton Size: 37x17x20cm



HS-M03 Dental Cotton Rolls

Material: 100% Cotton, natural. super Absorbent

Size(cm)	pcs/bundle	bundles/bag	bags/ctn	Carton Size(cm)
0.8x3.8	50	20	20	50x32x40
1.0x3.8	50	20	20	60x38x40
1.2x3.8	50	20	10	43x37x40
1.4x3.8	50	20	20	50x32x40



HS-M04 Dental Tips

Dia: 2.5mm, 2.0mm, 1.5mm 1.2mm
Packing: 100pcs/bucket, 4buckets/box, 40boxes/ctn
Carton Size: 52x38.5x228cm



HS-M05 Dental Needle

Parts: Stainless steel casing, AISI Standard 304
Medical grade PVC, within the spiral
Automatic out-needle device
Top surface: Opaque non-slip PE
Size: 25Gx25mm, 25Gx32mm, 25Gx38mm, 27Gx25mm, 27Gx32mm, 27Gx38mm, 30Gx25mm, 30Gx32mm, 30Gx38mm
Packing: 1pc/plastic bag, glue paper indicate size
50pcs/box, 100boxes/ctn 40x30x33cm

HS-M06 Dental Sponges

Material: Cotton
Cotton Yarn: 40'sx40's, 32'sx32's
Mesh: 19x9, 20x12 etc.
Ply: 8ply etc
Sterile or non-sterile

Size	pcs/pack	packs/ctn	Carton Size(cm)
2"x2"-8 ply	100	50	52x42x45
3"x3"-8 ply	100	50	52x50x30
4"x4"-8 ply	100	50	52x52x40



HS-M08 Dental Irrigation Needle

Material: Stainless Steel
Type: Bend needle, double vented, side vented, notched
Size: Bend needle: 18G, 20G, 21G, 22G, 23G, 24G, 25G
Length: 25mm
Sterile: EO
Packing: 100pc/PE Bag, 100bags/box, 6boxes/ctn
Carton Size: 68x32x38cm



HS-M07 Dental Tray

Content Of The Set:

Plastic tray,
Periodontal probe with plastic handle,
Plastic-handle tweezers or stainless metal points or a metal-handle
tweezers or forceps,
Plastic stomatoscope, glass mirror or pc mirror,
Abs handle, Dental bib,
Additional supplies such as a saliva ejector, gel tray, gauze swabs,
cotton balls, dental brush, dental cotton, mask, etc.
Can be assembled to the individual's specifications.

Sterilized with E.O. Gas.

Expiration Date: 2 Years

Packing: 1 Set/Plastic bag, 200sets/ctn 42x36x26cm



HS-M09 Dental Barrier Film


Material: PE

Core material: Plastic,paper

Color: Blue,transparent,customized

Size: 4"x6"(10.16cmx15.24cm) Thickness: 50micron

Packing: 1200pcs/roll, 1roll/box, 8box/ctn 31x30x26cm



HS-M11 Dental Elastomeric Power Chain

Material: Non-latex material

Type: short, long, closed

Color: Colorful

Packing: 1200rolls/ctn 46x30x31cm



HS-M10 Dental Ligature Tie

Material: Non-latex material, has continuous elastic

Type: long distance, short distance

Color: Colorful

Packing: 40pcs/bag,400bags/ctn/46x30x31cm



HS-M13 Dental Syringe

Material: Medical grade PVC

Color: White, transparent

Size: 5ml,12ml

Packing:

5ml 50pcs/bag, 2400pcs/ctn 42x42x36cm

12ml 50pcs/bag, 1250pcs/ctn 42x42x36cm

HS-M12 Denture Box

Material: PP

Size: 85x82x28mm

Packing: 1pc/OPP bag,450pcs/ctn/52x42x45cm



HS-M14 X-Ray Barrier Envelopes

Material: PVC

Size: 22x35, 20x40, 31x41, 27x54, 57x76mm, available

Color: Black and clear

Packing: 300 or 100 or 50pcs/box (according to the size you order)

HS-M15 Anti-Fog Face Splash Shield

Material: Sponge+PET shield

Size: 32x19cm.32x22cm,34x26cm

Thickness: 0.18mm

Packing: 1pc/bag, 200pcs/ctn 69x35x39cm



Dental Cotton Roll

100% pure cotton
Sterile/Non sterile
10 times' absorbency,sink time less than 10s
No lint,remain tight even get wet
Beautiful Printing

Size #1 0.8cm x 3.8cm (0.315" x 1-1/2")
Size #2 1.0cm x 3.8cm (0.375" x 1-1/2")
Size #3 1.2cm x 3.8cm (0.472" x 1-1/2")
Size #4 1.5cm x 3.8cm (0.551" x 1-1/2")
L 1-1/2" x 1/2" M 1-1/2" x 3/8" S 1-1/2" x 5/16"
8 x 25mm 10 x 30mm

*Other dimension, weight and package can be produced as per requirements



Cotton Filled(Dental) Sponge

100% pure cotton
USP VII gauze covered out side,cotton inside
Sterile/Non-sterile
Paper pouch/Paper-film pouch/Blister pack /Paper sleeve
Beautiful printing
High Absorbency and low lint

2" x 2"-8 ply 3" x 3"-8 ply 4" x 4"-8 ply

*Other dimensions or package can be produced as per requirements



Face Mask

PP + filter paper + pp
Paper/paper-film pouch or blister pack
Ear loop or tie on
White/Blue /Green /Pink /Yellow
Sterile/non-sterile

175mm x 95mm 145mm x 95mm

*Other weight, design and package can be produced as per requirements



Dental Bibs

Quality Level: A ; B; C
Paper: convention:16-17 gsm; customize:17-18gsm;18-19gsm
Film: covention:12-13 gsm
Pattern: convention: point; customize: can imprint the pattern and logo
Layer: convention: 2 ply paper+1 ply film; customize :1 ply paper+1 ply film



Tongue Depressor

Chinese white birch wood
Straight edge/round edge
Natural smell, taste free, clean, smooth
Sterile/Non-sterile
Bulk in box or individually packed in paper pouch

150*18*1.6mm 150*19*1.6mm 140m*14 *1.6mm
140*18*1.6mm 150*20*1.6mm

*Other dimensions or package can be produced as per requirements



Sterilization Pouch

Medical Grade Paper(60GSM or 70GSM)
Transparent Composite Film(White, Blue or Green)
Adopt Environmental Printing Ink
With STEAM, EO Chemical Indicator
Apply to hospital, laboratory, dental clinic nail,
beauty supply,tattoo&piercing supply etc.

50x200mm, 55x200mm, 75x200mm,100x200mm, 150x200mm, 200x200mm
250x200mm, 300x200mm, 350x200mm, 400x200mm



FOR ORDERS
Email: supplies@kumis.co.za
www.kumis.co.za