

BANDAGE/GAUZE





Tubular bandage (viscosa)

This kind of tubular bandage are soft .cool and elastic ,with a green line(or other color) in the middle of the bandage ,different size is available.

Tubular bandage (latex free)

Tubular bandage latex free ,with white color and skin color , different size is available .fit to feet. arm. elbow and so on.



Tubular bandage (latex)

This kind of tubular bandage have a good elasticity ,plain and soft ,different size is available fit to feet, arm, elbow and so on.

Stockinette

It is the best dressing product of tubular bandage ,high elasticity and good breathability .even in pressure ,tight to skin , different size is available.





100% cotton elastic bandage

100% cotton elastic bandage is the main product of our factory ,we are one of the largest manufacturers to produce the bandage in China

——materials: 100% high quality cotton

——color: natural White, bleached White, flesh, green, with or without blue/Red lines.

——weigh: 60–200g/m²

We can supply bandage with width ranging from 2.5cm to 30cm and length according to customers requirement .



100% cotton plain elastic bandage
(raw edge)

100% cotton plain elastic bandage have good elasticity ,good recovery ratio,good appearance and easily using. different size & color is available.

100% cotton plain elastic bandage
(woven edge)

100% cotton plain elastic bandage have good elasticity ,good recovery ratio,good appearance and easily using. different size & color is available.



100% cotton crepe elastic bandage
(raw edge)

100% cotton crepe elastic bandage have high absorbency and high elasticity, usually use in wound. surgical .clinical. medical usage ,washable . different size & color is available.



Spandex elastic bandage

Spandex elastic bandage is made of natural cotton and spandex fiber with fixed ends, comfortable . for medical clinic .health care and athletic spots use .

- Width within 2.5–20cm and the length according to customs requirements
- Weight: 60–200g/m²
- Packing: each roll with clips, individual pack then a dozen in a bag or box, several bag/box in a carton
- Color: Nature White, bleached White flesh and so on



Spandex elastic bandage

spandex elastic bandage are made by cotton yarn and spandex threads , friendly to people ,to fix wound dressings for first aid. High elasticity ,washable ,sterilizable .



Spandex elastic bandage (plain & crepe)

spandex elastic bandage are made by cotton yarn and spandex threads , friendly to people ,to fix wound dressings for first aid. High elasticity ,washable ,sterilizable .

High elastic bandage

High elastic bandage is cool and comfortable to wear ,superior strength and elasticity ,with selvages and fixed ends.

- materials: cotton.polyester.rubber.
 - packing :each bandage is wrapped individuauly in waterproof bag,the outer package is strong cardboard carton to keep best storing condition .
 - weight : 60–200g/m²
 - different width and length available
- For more information, pls contract us





High elastic bandage(latex)

Latex elastic bandage are orthopedic grade woven bandages that provide consistent elasticity, support and compression, made of the finest quality yarns, these bandage can withstand numerous washing and sterilization, the material of latex rubber is gain authentication, that is innocuous to the skin of patient.



High elastic bandage(latex free)

Latex free high elastic bandage, orthopedic grade knitted bandage that can be wrapped and secured easily, it provides patient comfort and is durable enough to be reusable. For use in hospitals, physician office, school health and home care.



PBT bandage

PBT bandage (conforming gauze bandage) , it applies to the orthopaedic surgeons, it feels comfortable 、lightness and elegant. good convenient, excellent elastic and reasonable price.

——It made of polyester fiber 、 polyamide 、 cotton or DTY

——Width within 2.5~20cm and the length according to customs requirements

packing: individual packing, then a dozen in a bag /box



PBT elastic bandage (thin)

PBT elastic bandage , plain or crepe . used to treat wound care or the daily accident .It's made of soft texture ,with high and permanent elasticity



PBT elastic bandage (thick)

PBT elastic bandage (thick) made of soft texture and cotton , with high and permanent elasticity ,suitable for the use of binding all parts of body, with woven edges, with smooth surface,200% elasticity ,it is ideal for fixation of wound dressing.



Elastic bandage

This kind of bandage made of of soft texture ,spandex and cotton. Soft , elasticity 280~290% ,high absorbency. washable ,sterilizable,color can be white and flesh etc .

Gauze Bandage

Gauze bandage (w o w bandage) cutting edges or woven edges , it has high absorbent and soft ,pure white ,widely used and cut out easily ,we can offer all kinds of gauze products with good quality and competitive price .

——material :pure absorbent cotton 21's .32's .40's

——mesh: 13 \17\20\24 threads





Gauze bandage (woven edges)

side woven gauze bandage is made of 100% cotton with folded edges, soft and comfortable ,more convenient and more advantage.we can supply gauze bandage width within 1cm–20cm,and any length according to customers .



Gauze bandage (cut edges)

this gauze roll are made of pure cotton yarn ,through high temperature and pressure degreased and bleached ready–cut, it is the necessary products for hospital, White paper per roll .blue kraft per dozen .



First aid bandage

First aid bandage usually used for emergency aid , it consists of absorbent pad and elastic bandage, soft and comfortable to the wound, easy to operate, various size for the bandage and the dressings for your choice.



First aid bandage(thick)

First aid bandage unite the elastic bandage and the wound dressing, integrating the function of wrapping , stoping bleeding and accessorial fixing ,it can be used in different injures in different parts of body in case of battlefield or pre-hospital emergency.



Plaster of paris bandage

Plaster of paris bandage is suitable for instant use , plaster covers the gauze completely, it is necessary product for surgery treatment.
Width:5cm/7.5cm/10cm/12.5cm/15cm/20cm, length: 2.7m/3.0m/4.6m



Cohesive elastic bandage

Cohesive elastic bandage , excellent pressure bandage , without cutting off circulation , self adhesive to itself , no clips or fasteners needed , easily hand tearable . its made of nonwoven or cotton .different color & size is available .



First aid bandage(thin)

Usually made of PBT bandage & absorbent pad. There is a thin membrane-like material in non-ticks pad on the outer side toward the wound ,it allows the body fluid to fully pass to the absorbent layer , it also makes it easier to take the dressing from the wound . Available in multiple package , non sterile and sterile , individual pack in cellophane or pouch .



Orthopaedic padding

Orthopaedic padding is wide-used supporting material for surgery and wadding of plaster of paris bandages, it is soft ,smooth and comfortable to body's shape .its made of cotton and polyester,
Width:5cm/7.5cm/10cm/12.5cm/15cm/20cm, length:2.7m/3.0m/4.6m



Kerlix Bandage

100% pure cotton
High absorbency, great softness
Sterile or non sterile

4.5" x 4.1 yards-6 ply

*Other dimensions and type can be produced as per requirements



First Aid Compress Bandage

Weave/knitted, plain/crepe
Sterile or non sterile
Individually packed in cellophane, pouch or blister
Different weight and material

Stretched Length:
4cm x 4m 6cm x 4m 8cm x 4m
10cm x 4m 12cm x 4m 15cm x 4m

*Other dimensions, weight or package can be produced as per requirements



High Elastic Bandage

Weave/knitted, bleached/un bleached, skin color
Individually packed in cellophane
Different weight
Cotton/spandex

Stretched Length:
2.5cm x 4m 5cm x 4.5m 7.5cm x 4.5m
10cm x 4.5m 15cm x 4.5m

*Other dimensions, weight or package can be produced as per requirements



Crepe Bandage

Weave/knitted, bleached/un bleached, skin color
Sterile or non sterile
Individually packed in cellophane, pouch or blister
With/without red/blue thread
Different weight
Cotton/spandex

Stretched Length:
2.5cm x 4m 5cm x 4.5m 7.5cm x 4.5m
10cm x 4.5m 15cm x 4.5m

*Other width, length, weight and package can be produced as per requirements



Stockinette

Natural white or bleached white
100% cotton
Ideal for sprains, strains, joint support and dressing retention

5/7.5/10/15/20/25cm x 5/10/25m

Packing: each piece in one polybag or despende box

*Other width, length and package can be produced as per requirement



Plaster of Paris Bandage

Consists of plain-weave gauze fabric coated with blend of alpha and beta calcium sulphate crystals, spooled on round perforated plastic core

Available in different size
Available in white or colored
Immersion time 2 to 3 seconds only
Excellent moulding ability

*Other specifications or package can be produced as per requirements



PBT Conforming Bandage

Weave/knitted, plain/crepe
Sterile or non sterile
Individually packed in cellophane, pouch or blister
Different weight and material

Stretched Length:
4cm x 4m 6cm x 4m 8cm x 4m
10cm x 4m 12cm x 4m 15cm x 4m
1" x 4.1 yards 2" x 4.1 yards 3" x 4.1 yards
4" x 4.1 yards 6" x 4.1 yards

*Other width, length, weight or package can be produced as per requirements



Triangular Bandage

BP 13, 17, 20 threads gauze or Non-woven
Bleached or unbleached
Individually packed with 2 safety pins

37" x 37" x 52" 40" x 40" x 56"

*Other dimension, weight or package can be produced as per requirements



Gauze Bandage

BP 13, 17, 20 threads, white/green or blue gauze
Sterile or non sterile
12pcs/paper pack or individually packed
Cutting edges/Selvage

1"/2"/3"/4"/6" x 4/10 yards

*Other width, length and package can be produced as per requirements



Absorbent Gauze Roll

USP VII,BP 13,17,20 threads , white/green or blue gauze
With/without X-ray detectable thread
1-ply,2-ply or 4-ply

Jumbo roll 36"/48" x 1000/2000 m
36"/48" x 5.5/11/55/100/110 yards

*Other dimensions,width,length and package can be produced as per requirements



Gauze Swabs

100% pure cotton USP VII,BP 13,17,20 threads , white,green or blue gauze
With/without folded edge
With/without X-ray detectable thread
Available in STEAM,EO OR GAMMA sterile

5 x 5cm-8/12/16 ply
7.5cm x 7.5cm-8/12/16 ply
10cm x 10cm-8/12/16 ply
10cm x 20cm-8/12/16 ply

100/200 pcs/paper sleeve or box
1/2/3/4/5/10 pcs/pouch or blister

*Other dimensions,plies, and package can be produced as per requirements



Laparotomy Sponge

BP 13,17,20 threads , white/green or blue gauze or non-woven
With/without X-ray detectable thread or chip
With/without metal or plastic ring,blue or white cotton loop
Pre-washed or unwashed
100pcs/polybag
5pcs/blister with/without white or blue medical wrapping paper

20cm x 30cm-4/6 ply	30cm x 30cm-4/6 ply
45cm x 45cm-4/6 ply	50cm x 60cm-4/6 ply
90cm x 8cm-4ply	40cm x 40cm-4 ply
20cm x 20cm-12 ply	30cm x 30cm-12 ply
45cm x 45cm-12 ply	
18" x 18"-4ply	12" x 12"-4ply
4" x 18"-4ply	8" x 36"-4ply

*Other dimensions,plies and package can be produced as per requirements



Y/I Cut Gauze Swabs

100% pure cotton USP VII,BP 13,17,20 threads

Y-cut or I-cut

Sterile or non sterile

100/200 pcs/sleeve or in paper-paper/paper-film/blister pack


100% cotton swab:

10cm x 10cm-12/16 ply

Non-woven swab:

5cm/10cm-6 ply

*Other dimensions ,plies and package can be produced as per requirements



Gauze Cutting

100% pure cotton USP VII,BP 13,17,20 threads , white,green or blue gauze

With/without X-ray detectable thread

Sterile/Non sterile

Packed in polybag/paper-paper/paper-film or blister pack

10cm x 10cm 10cm x 20cm 20cm x 20cm

25cm x 25cm 30cm x 30cm 30cm x 40cm

*Other dimensions can be produced as per requirements



Gauze Ball

100% pure cotton BP 17,20 threads

With/without X-ray detectable thread

In round,flat or hat shape

Sterile or non sterile

Unfolded size:

10cm x 10cm 15cm x 16cm 20cm x 20cm

25cm x 26cm 30cm x 30cm 30cm x 40cm

40cm x 50cm

*Other dimensions and package can be produced as per requirements



Gauze Sheet

100% pure cotton USP VII,BP 13,17,20 threads , white,green or blue gauze

Sterile/Non sterile

Packed in polybag/paper/paper-film bag or blister pack

90cm x 18m 90cm x 27m 90cm x 36m

90cm x 1m 90cm x 3m 90cm x 5m

90cm x 10m

1 square meter 1/2 square meter

1/4 square meter

*Other dimension or package can be produced as per requirements



Post-op Sponge

Covered with gauze/non woven
Very high absorbency while with no lint
Tissue in the middle, cotton inside or tissue inside only
Machine/Manual folding
Sterile/Non sterile
Paper/paper-film pouch /blister pack

3" x 3" 3" x 4" 4" x 4"

*Other dimension and package can be produced as per requirements



Vaseline Gauze

Consists of gauze, paraffin/vaseline, oil emulsion
21S/24X24
200±20g/m²
Available in different size
1pc/individual aluminum foil bag packing; 5pcs to 100pcs/mid box

5x5cm, 10x10cm, 10x20cm, 10x40cm,
10x700cm, 15x100cm, 15x200cm

*Other specification can be produced as per requirements



Non-woven Swabs

70% viscose + 30% polyester
30gsm, 35gsm, 40gsm, 45gsm or 50gsm
With/without X-ray detectable thread
High absorbency and low lint

5 x 5cm-4/6/8 ply
7.5cm x 7.5cm-4/6/8 ply
10cm x 10cm-4/6/8 ply

1/2/3/4/5/10 pcs/pouch or blister
100/200 pcs/sleeve

*Other package can be produced as per requirements



Drain Sponge

100% pure cotton USP VII, BP 13, 17, 20 threads
or 30/35/40/45/50 gsm non-woven
Y-cut or I-cut
Sterile or non sterile
100/200 pcs/sleeve or in paper-paper/paper-film/blister pack

100% cotton swab:
10cm x 10cm-12/16 ply
Non-woven swab:
5cm/10cm-6 ply

*Other dimensions ,plies and package can be produced as per requirements



Non-woven Ball

100% pure cotton , 30/35/40/45/50 gsm non-woven
With/without X-ray detectable thread
In round, flat or hat shape
Sterile or non sterile

Unfolded size:

10cm x 10cm	15cm x 16cm	20cm x 20cm
25cm x 26cm	30cm x 30cm	30cm x 40cm
40cm x 50cm		

*Other dimensions and package can be produced as per requirements



Non-woven Lap Sponges

30gsm, 35gsm, 40gsm, 45gsm or 50gsm
With/without X-ray detectable thread or chip
With/without metal or plastic ring, blue or white cotton loop
Pre-washed or unwashed
100pcs/polybag
5pcs/blister with/without white or blue medical wrapping paper

20cm x 30cm-4/6 ply	30cm x 30cm-4/6 ply
45cm x 45cm-4/6 ply	50cm x 60cm-4/6 ply
90cm x 8cm-4ply	40cm x 40cm-4 ply
20cm x 20cm-12 ply	30cm x 30cm-12 ply
45cm x 45cm-12 ply	
18" x 18"-4ply	12" x 12"-4ply
4" x 18"-4ply	8" x 36"-4ply

*Other dimensions, plies and package can be produced as per requirements




Abdominal Pad(ABD)

Non woven cover outside, tissue and cellulose in side
All edges are sealed to prevent linting
Sterile/Non-sterile

5" x 9" 7.5" x 8" 8" x 10"

*Other dimensions or package can be produced as per requirements



Underpads

Top layer: non-woven fabric and bottom is PE film
Middle layer: consists of imported fluff pulp and SAP
Rapid absorbency can absorb the pets' urine,
keep house and garden clean and dry, and prevent from bad smell

60x40cm, 60x45cm, 60x60cm,
75x45cm, 75x60cm, 90x60cm 80x150cm etc

6040, 6045, 6060, 7545, 7560, 9060cm
According to customers requirements, L:45to180cm, W:30to80cm

10pcs/bag, plastic internal packing; 10 bags/carton



Absorbent Cotton Wool

100% pure cotton, high absorbency
50g/100g/200g/250g/400g/454g/500g/1000g/50kg
Wrapped by blue/white medical paper or in polybag
50kg in big bale, other sizes in cartons

*Other weight and package can be produced as per requirements



Cotton Ball

100% pure cotton, high absorbency
Sterile/Non-sterile
Packed bulk in polybag, sterile product in paper bag/blister pack

0.2g/0.3g/0.4g/0.8g/1.0g/2.0g

*Other weight and package can be produced as per requirements



Cotton Pad/Wipe

100% pure cotton, great softness and high absorbency
Different size/gramme weight
Sterile/Non-sterile
Paper/paper-film pouch /blister pack

4cm x 4cm 5cm x 5cm
7.5cm x 7.5cm 10cm x 10cm

*Other dimension, weight or package can be produced as per requirements



Cotton Tipped Applicator

100% cotton, high absorbency
Wooden stick/plastic stick
Single tip/Dual tip, different sizes
Sterile/Non-sterile

3"/6" x 5/7 mm

*Other specification and package can be produced as per requirements



Undercast Padding/Orthopedic Bandage

100% pure cotton/Synthetic
Sterile/Non-sterile
Paper pouch/Paper-film pouch/Blister pack
Different width and length

2"/3"/4"/6" x 4 yards

*Other width ,length or package can be produced as per requirements



Cosmetic Cotton

100% pure cotton,great softness and high absorbency,no lint,no irritation
0.4g/0.55g/0.6g
In square type or round type
50/80/100/160 pcs Packed in polybag or box
Excellent for cleaning make-up,removing lipstick,removing nail polish,
applying moisturizers,patting with cleansing lotion and so on

*Other dimension ,weight or package can be produced as per requirements



X-ray Peanuts

100% pure cotton BP 17,20 threads
X-ray Detectable
In Peanuts Shape
Sterile or non sterile

Unfolded size:

10cm x 10cm	15cm x 16cm	20cm x 20cm
25cm x 26cm	30cm x 30cm	30cm x 40cm
40cm x 50cm		

*Other dimensions and package can be produced as per requirements



Tubular bandage

Tubular bandage (stockinette bandage) is made of 100% cotton or with spandex , viscose other elastic fibers, soft and comfortable, different size is available.fit to feet.arm.elbow and so on.

Size,width(natural) application

A	4.50cm	infant feet or arms
B	6.25cm	small hands and arms
C	6.75cm	medium arms small ankles
D	7.50cm	large arms .medium ankles.small kness
E	8.75cm	large ankles.medium knees.samll thighs
F	10.00cm	large knees、 medium thighs
G	12.00cm	large thighs
J	17.50cm	small trunks
K	21.50cm	medium trunks
L	32.50cm	large trunks
M	37.50cm	extra-large trunks

We can supply specified size from clients.

Width : 1.5cm~33cm any length.

Packing: each roll in a OPP bag or box with print.



FOR ORDERS
Email: supplies@kumis.co.za
www.kumis.co.za