

HOSPITAL BEDS/ CABINETS/TROLLEYS



HOSPITAL BEDS

ICU – WARD - GENERAL USE - HOME CARE

ICU - MULTIFUNCTION ELECTRIC BED



Standard Features: ICU Bed (Five function) with electrically operated remote controlled back rest tilting 0-70°, knee rest tilting 0-35°, Trendelenburg tilting 0-12°, Reverse Trendelenburg tilting 0-12° Size: L 2160mm x W 970mm x H 520-920mm (Adjustable Height)

* ELECTRIC MODELS AVAILABLE *

Detailed information will be supplied on request

- Multifunction ICU Bed- Electric
- Five Function – Electric
- Three Function- Electric

HOSPITAL BEDS/CABINETS/TROLLEYS

HOSPITAL BEDS

ICU - WARD - GENERAL USE - HOME CARE

Standard Care Bed

Platform base -Mesh/Plate

Head and Foot bows

Accessories:

Safety sides

Castor Wheels

Mattress



Single Function Care Bed/Pneumonia Bed

2 Section Platform base -Mesh/Plate

Head and Foot bows

Manual Backrest

Accessories:

Safety sides

Castor Wheels

Mattress



* FEATURED ITEM: ISOLATION BED



Simple and basic design isolation bed

Optional:

Safety sides

Mattress

Castors

HOSPITAL BEDS

ICU – WARD - GENERAL USE - HOME CARE

"Hospital beds are one of the most effective types of beds for anyone with mobility issues. And, they are not just for hospitals"



SINGLE FUNCTION BED

- *Detachable ABS head/foot board, steel bed plate, powder coated
- *Back lifting: 0 ~ 80° (±5°)
- *Aluminum alloy side rails

TWO FUNCTION BED

- *Detachable ABS head/foot board, steel bed plate, powder coated
- *Back lifting: 0 ~ 80° (±5°)
- *Knee lifting: 0 ~ 40° (±5°)
- *Aluminum alloy side rails



THREE & FIVE FUNCTION BED

- *Aluminum alloy side rails
- *Back lifting: 0 ~ 80° (±5°)
- *Knee lifting: 0 ~ 40° (±5°)
- *Height adjustment: 420~680mm
- *Trendelenburg & Rev.trendelenburg: 0~12°
- * 5" universal durable castors, two with brake
- * With I.V.rod

HOSPITAL BEDS/CABINETS/TROLLEYS

PEDIATRIC AND NEONATAL CARE

INFANT & CHILD COTBEDS



"Pediatric beds are designed for infants, toddlers, and young children
Also known as a safety bed or enclosed bed."

FEATURES AVAILABLE

Backrest – Gas bottle holder
Trendelenburg Function – IV Pole

BASSINETTES

A healthcare bed designed for infants, with a removable tub- **Stainless steel optional**



Adjustable by use of a gas lift



Fixed Model /Storage

ACCESSORIES FOR CARE

Gas bottle holder- Iv Pole – Transparent crib- Mattress- lower storage Cabinet

SWADDLING TABLE

Upholstered top with steel frame on wheels




CSSD - LOGISTICS

Shelving – Packing tables – Sorting tables -Trolleys

Packing tables Transport trolleys Storage solutions



The central sterile services department (CSSD), is an integrated place in hospitals and other health care facilities that performs sterilization and other actions on medical devices, equipment and consumables

ITEM	PICTURE	PRODUCT CODE	PRODUCT NAME
1		PTSS 1400 X 600-900 x900H PTSS 1600 x 600-900 x900H PTSS 1800 x 600-900 x900H PTSS 2000 x 600-900 x900H	Stainless steel Packing Table
2		TT 720 Trolley for wrapping paper	Stainless steel wrapping paper and transportation trolley
3		OTT -S * Single Open trolley OTT -D * Double open trolley	Stainless steel Open Transport trolley Compatible with wire baskets or shelves
4		SS5F - Fixed - Made with solid shelves stainless steel SS5O -Open - Made with tubular supports for ventilation	Stainless steel tubular shelving with solid shelves or tubular shelving.
5		SS5W -Wire - Made with Wire shelves	Stainless steel tubular shelving made with wire shelves

HOSPITAL BEDS/CABINETS/TROLLEYS

CSSD - LOGISTICS

Shelving – Packing tables – Sorting tables -Trolleys

STERILE AND NON-STERILE TRANSPORTATION TROLLEY



Made from stainless steel
2 x shelves with ventilation
Mounted on 4 x 125mm castors



HEAVY DUTY INSTRUMENT TROLLEY FOR CSSD

Heavy-duty stainless-steel trolley for
Bulk transport on instruments
Bi-layer design 2 & 3 Shelves Variation
Standard Size: 900mm x 600mm
4 x 100mm Castors

PACKING AND SORTING TABLES

Made from stainless steel
Bi-layer design
Standard size 1800x900x900
Optional
4 x 100mm castors



LAUNDRY - LOGISTICS

LAUNDRY LOGISTICS



Soiled linen receiver with bag
Available in epoxy coated & Stainless steel
Models Available: **Single Bag trolley**
: **Double Bag trolley**



Clean & Soiled linen trolley with 3 shelves & bag
Available in epoxy coated & Stainless steel



Clean Linen Trolley with 3 shelves
Available in epoxy coated & Stainless steel

Soaking Bin Trolley
Bin colors: Red-Blue-White-Red



Towable Bin trolley on galvanized frame
Bin colors: Red-Blue-White-Red

HOLLOWWARE



- 15cm Kidney Dish, Stainless steel
- 20cm Kidney Dish, Stainless steel
- 25cm Kidney Dish, Stainless steel
- 30cm Kidney Dish, Stainless steel
- 15cm Kidney Dish, Plastic
- 20cm Kidney Dish, Plastic
- 25cm Kidney Dish, Plastic
- 30cm Kidney Dish, Plastic
- 50ml Gallipot, Stainless steel
- 90ml Gallipot, Stainless steel
- 120ml Gallipot, Stainless steel
- 180ml Gallipot, Stainless steel
- Bedpan, stainless steel
- Bedpan, Fracture type stainless steel
- Bedpan, Fracture type Plastic
- Urinal, 1 Liter, stainless steel, jug type
- 1 Liter Measuring Jug, Stainless steel

Our Hospital Hollowware range provide great convenience for multiple uses during medical procedures, surgery & patient care.



* Plastic & Stainless steel available

THEATRE AND SLUICE ROOM EQUIPMENT

BOWL STANDS & KICKABOUTS

- Single Bowl stand
- Double bowl stand
- Kickabout – Bowl
- Kickabout - Bucket



Made from stainless steel

STERILE - INSTRUMENT STORAGE AND TRANSPORT TROLLEY



Made from Stainless Steel
Non-Marking Rubber Castors
3-5 Wire shelves

Made from stainless steel

MAYO TROLLEY

Made from Stainless Steel
Non-Marking Rubber Castors
Fixed & hydraulic



THEATRE TROLLEYS

INSTRUMENT TROLLEY

Stainless steel Sizes

457 mm (w)

610 mm (w)

915 mm (w)



OPTIONAL: FIXED shelves

DRESSING TROLLEY

Stainless steel Sizes

457 mm (w)

610 mm (w)

915 mm (w)



OPTIONAL: FIXED SHELVES

DRESSING TROLLEY WITH BUCKET AND BOWL



Stainless steel
Size 915 mm (W)

MOBILE IV STANDS

5 LEGGED MOBILE IV STAND

- **Stainless Steel & Epoxy coated options available**
- Nylon or chrome base options based on the environment
- 2IV ,4IV & 6 IV Hooks

Weighted base options available



NESTING MOBILE IV POLES

Stainless Steel & Epoxy coated options available
Used for easy storage



*More iv stands under the general ward furniture category

INTENSIVE CARE TROLLEYS

CT1– Chart trolley used for neonatal and infant ICU

1 drawer and 1 basket with runners.

- Optional is 2 drawers



BEU-1- ICU bed end unit lockable. Fully enclosed with chart holder on 1 side

- Wooden top in choice of colors.
- Optional tops: stainless steel, hard plastic.



CT-2- Chart trolley with two drawers. Simple and economical.

- Open frame trolley. Lockable draws.
- Optional tops: stainless steel, hard plastic



BEU-2- ICU bed end unit lockable.

Open frame with chart holder on 1 side. Wooden top in choice of colors. Drawers include runners.

- Optional tops: stainless steel, hard plastic.



PATIENT TRANSPORTATION

PATIENT EMERGENCY STRETCHER

Product Functions:

- 1) Backrest max upward angle: 75°
- 2) Footrest Max upward Angle: 45°
- 3) Height adjustment: 720-950mm
- 4) Trendelenburg: 15°
- 5) Anti-Trendelenburg: 15°



PATIENT TRANSPORT TROLLEY - HYDRAULIC

Product features

- Hydraulic foot pedal controlled Hi-Lo position
- Tilting assisted by gas spring
- Backrest adjustable by gas spring
- PVC bumper guard rail
- Push handle at head and foot ends
- Collapsible side rails
- IV pole
- Storage utility tray
- Oxygen tank holder



PATIENT TROLLEY - FIXED HEIGHT

Product features

- Adjustable backrest by gas spring
- Safety sides
- Push handle at head and foot ends
- Collapsible side rails
- IV pole
- Storage utility tray
- Oxygen tank holder



CADAVER TROLLEYS



MORTUARY TROLLEY - WITH CONCEALMENT

Epoxy coated frame with stainless steel body tray
Available in stainless steel

CADAVER TROLLEY WITH REMOVABLE BODY TRAY

Fully Stainless steel
Available with concealment



HYDRAULIC LIFTING TROLLEY - STAINLESS STEEL

Fully stainless-steel construction
Powered by manual hand pump
Removable body tray



MORTUARY/AUTOPSY TABLE

Fully stainless-steel construction 2200x800x900

4 x 125mm castors

MEDICINE & ANESTHESIA

MEDICINE TROLLEY - ABS

Made from ABS and aluminum

4 & 6 Drawer combinations

**Also available in steel with epoxy coating*

Detailed information on request

"mainly consists of groups of drawers for storing different kinds of medicines"



ANESTHESIA CART - ABS

With Overhead Storage Bins

Made from ABS and aluminum

4,5 & 6 Drawer combinations

**Also available in steel with epoxy coating*

Detailed information on request



"Anesthesia carts are hospital devices used to store tools that are necessary for aid during procedures that require administration of anesthesia"



FOR ORDERS
Email: supplies@kumis.co.za
www.kumis.co.za