

INFUSION PRODUCT





HS-D01 Disposable Infusion Set

Material: Medical grade PVC
Components: Small Chamber+plastic needle+flow regulator+ rubber gum bubble+with Latex Tube+with Filter
Color: Transparent
Sterilization: EO
Packing: 1pc/PE Bag or blister bag, 500pcs/ctn
Carton Size: 58x38x38cm



HS-D101 Disposable Infusion Set

Material: Medical grade PVC
Components: Small Chamber+stainless needle+flow regulator +rubber gum bubble+with Latex Tube+with Filter
Color: Transparent
Sterilization: EO
Packing: 1pc/PE Bag or blister bag, 500pcs/ctn
Carton Size: 58x38x38cm



HS-D02 Disposable Infusion Set

Material: Medical grade PVC
Components: Small Chamber+plastic needle+flow regulator+with Y site+with Filter+with luer lock adapter
Color: Transparent
Sterilization: EO
Packing: 1pc/PE Bag or blister bag, 500pcs/ctn
Carton Size: 58x38x38cm



HS-D03 Disposable Infusion Set (Dark Color)

Material: Medical grade PVC
Components: Small Chamber+plastic needle+flow regulator +rubber gum bubble+ with Latex Tube+with Filter
Color: Brown
Sterilization: EO
Packing: 1pc/PE Bag or blister bag, 500pcs/ctn
Carton Size: 58x38x38cm



HS-D08 Disposable Infusion Set With Regulator


Material: Medical grade PVC
Components: Small Chamber+plastic needle+flow regulator+with Y site+IV regulator+with Latex Tube+with Filter
Color: Transparent
Sterilization: EO
Packing: 1pc/PE Bag or blister bag, 300pcs/ctn
Carton Size: 61x38x28cm



HS-D05 Burette Infusion Set

Material: Medical grade PVC
Components: Air filter+soft tube+drip chamber+flow regulator+fluid filter+clamp+stop valve+Y injection site +protective cap
Size: 100ml, 150ml, 20 drops/ml, 60drops/ml
Color: Transparent
Sterilization: EO
Packing: 1pc/PE Bag or blister bag, 100pcs/ctn
Carton Size: 65x38x33.5

INFUSION PRODUCT



HS-D06 Disposable Blood Transfusion Set

Material: Medical grade PVC

Components: Vented spike+drip chamber+membrane filter+flow regulator+latex tube+luer lock connector+needle

Options available:

With or without air vented spike,

With or without needle,

With or without "Y" injection port,

Luer lock or luer slip connector.

Color: Transparent

Packing: 1pc/PE Bag or blister bag, 500pcs/ctn

Carton Size: 58x38x38cm



HS-D106 Disposable Blood Transfusion Set With Regulator

Material: Medical grade PVC

Components: Vented spike+drip chamber+membrane filter+IV regulator+latex tube+luer lock connector

Options available:

With or without air vented spike; With or without needle,

With or without "Y" injection port,

Luer lock or luer slip connector.

Color: Transparent

Packing: 1pc/PE Bag or blister bag, 500pcs/ctn

Carton Size: 58x38x38cm



HS-D07 Stopcock Three Way

Material: PC+PE+ABS

Components: Body+plug+screw cap

Options available: with tube or without tube

Test pressure: 5.5kg.

Products Thread: 6:100 screw thread standard

This in line with international GB standards

Packing: 1pc/hard blister, 100pcs/box, 1000pcs/ctn

Carton Size: 47x28x54cm



HS-D10 Tubercle Bacillus Syringe

Material: Medical grade PP

Components: Barrel+Piston+Plunger+Needle (or without needle)

Type: Luer clip, luer lock

Size: 1ml

Packing: 1pc/PE bag or blister bag, 3600pcs/ctn

Carton Size: 16.5x48.5x38cm



HS-D09 Safety Auto Destruct Syringe

Material: Medical grade PP

Components: Barrel+Piston+Plunger+Clip+Needle (or without needle)

Size: 2ml, 3ml, 5ml, 10ml

Color: Transparent

Packing: 1pc/PE bag or blister bag

Carton Size:

2ml/3ml 3000pcs/ctn 59x44x45cm

5ml 2400pcs/ctn 56x44x46.5cm

10ml 1200pcs/ctn 56x44x46cm



HS-D11 Disposable Syringe

Material: Medical grade PP

Components:

2 parts: Barrel+Plunger+Needle (or without needle)

3 parts: Barrel+Piston+Plunger+Needle (or without needle)

Type: Luer slip, luer lock

Size: 1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 30ml, 50ml, 60ml

Packing: 1pc/PE bag or blister bag

Carton Size: PE bag packing as below:

1ml 3600pcs/ctn 58x45x44.5cm

2ml/3ml 2700pcs/ctn 58x45x44.5cm

5ml 2400pcs/ctn 58x45x44.5cm

10ml 1600pcs/ctn 58x43x44.5cm

20ml 900pcs/ctn 58x43x44.5cm

30ml 640pcs/ctn 58x43x44.5cm

50ml/60ml 300pcs/ctn 51x43x46cm



HS-D109 Irrigation syringe

Material: Medical grade PP
Components: Barrel+Piston+Plunger
Type: Cap or without cap
Size: 50ml, 60ml, 100ml, 150ml, 200ml
Sterilization: EO
Packing: 1pc/PE bag or blister bag
Carton Size: 50/60ml: 300pcs/ctn 51x43x46cm




HS-D209 Comestic Syringe

Material: Medical grade PP
Components: Barrel+Piston+Plunger
Color: Transparent
Size: 1ml, 3ml
Sterilization: EO
Packing: 1pc/PE bag or blister bag
Carton Size:
 1ml 3600pcs/ctn 58x45x44.5cm
 3ml 2700pcs/ctn 58x45x44.5cm



HS-D12 Insulin Syringe

Material: Medical grade PP
Components: Barrel+Piston+Plunger+Needle+Cap
Color: Orange
Type: 40IU, 100IU
Size: 0.3ml, 0.5ml, 1ml
Needle Size: 29G, 30G, 31G
Packing: 1pc/PE bag or blister bag
Carton Size:
 0.3ml/0.5ml 200pcs/box, 3000pcs/ctn 48.5x40.5x43.5cm
 1ml 100pcs/box, 3200pcs/ctn 61.5x48.5x38cm



HS-D13 Disposable Needle

Material: Stainless steel
Components: PVC base+Needle
Size: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G, 28G, 29G, 30G
Color: Pink, Cream, Yellow, Green, Black, Blue, Purple, Orange, Brown, Gray, Blue-green, Red, Light-yellow
Packing: 1pc/blister bag, 100pcs/box, 100box/ctn
Carton Size: 45x46x40cm



HS-D15 Spinal Needle

Material: Stainless steel
Components: PP hub+Needle
Type: Pencil point tip, Quincke tip
Size: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G
Color: Pink, Cream, Yellow, Green, Black, Blue, Purple, Orange, Brown, Gray
Length: 90mm-220mm
Packing: 1pc/blister bag, 1000pcs/ctn
Carton Size: 42x27x52cm



HS-D16 Epidural Needle

Material: Stainless steel
Components: PP hub+Needle
Size: 15-18G
Length: 80mm, 90 mm
Packing: 1pc/blister bag, 50pcs/box, 1000pcs/ctn
Carton Size: 54x42x36cm

INFUSION PRODUCT



HS-D17 Scalp Vein Set

Material: Medical grade PVC+Stainless steel
 Wing Type: Butterfly wing, Single wing
 Size: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G
 Color: Pink, Cream, Yellow, Green, Black, Blue, Purple, Orange, Brown, Gray
 Length: 15-30cm
 Packing: 1pc/PE bag or blister bag, 100pcs/box, 2000pcs/ctn
 Carton Size: 53x37x27cm



HS-D110 I.V. Cannula

Material: Stainless steel
 Type: With butterfly wings, Injection port with wings, Pen-like, Safety Pen-like
 Size: 14G, 16G, 18G, 20G, 22G, 24G, 26G
 Packing: 1pc/blister bag,
 Carton Size:
 Injection port: 1000pcs/ctn 72x31x31.5cm
 Pen-like: 1000pcs/ctn 47.5x38x25cm



HS-D18 Disposable Extension Set

Material: Medical grade PVC
 Tube Length: 10-25cm, 30-50cm, 80cm, 100-120cm
 Packing: 1pc/blister bag,
 Carton Size:
 10-25cm 50pcs/ctn 72x56x32cm
 30-50cm 50pcs/ctn 72x56x32cm
 80cm 50pcs/ctn 72x56x32cm
 100-120cm 40pcs/ctn 72x56x32cm



HS-D111 Safety I.V. cannula

Size: 18G, 20G, 22G, 24G, 26G
 Color: Green, Pink, Blue, Yellow, Purple
 Packing: 1pc/hard blister, 200pcs/box, 1200pcs/ctn
 Carton Size: 58x42x47cm



HS-D19 A.V. Needle


Material: Stainless steel
 Size: 15G, 16G, 17G
 Length: 30cm
 Packing: 1pc/PE bag or blister bag, 50pcs/box, 500pcs/ctn
 Carton Size: 63x36x45.5cm



HS-D309 Blood Bag

Material: Medical grade PVC
 Color: White transparent
 Packing: 1pc/PE bag

Specification Bag	Capacity	Anticoagulant	QTY	Carton Size
Single	250ml-500ml	CPDA/CPD	100	51x32x26cm
Double	250ml-500ml	CPDA/CPD	80	51x32x26cm
Triple	250ml-500ml	CPDA/CPD+SAG-M	30	36x36x20cm
Quadruple	250ml-500ml	CPDA/CPD+SAG-M	30	36x36x20cm



P1015001
Hypodermic Needle
Size : 15-31G



P1006001-LL
Disposable Syringe
Size:1/2/3/5/10/20/30/50/60ml



P1006001-LS
Disposable Syringe
Size:1/2/3/5/10/20/30/50/60m



P1015002
Disposable Syringe
Size:1/2/3/5/10/20/30/50/60ml



P1017001
Insulin Syringe
Size:0.3/0.5/1ml with needle 29-30G



P1018001
Disposable Syringe
SIZE: 29~32G



P1014001
Holder with Safty cover
Modle :standard



P1018001
Enfit Syringes
Size: 3/5/10/20/30/50/60/100ml



P1006003
Glass Prefilled Syringe
Size:1ml/3ml/5ml,with needle or luerlock



P1019001
Angiographic Syringe
Size :60/100/200ml



P1020001
Irrigation Syringe
Size:60/100/150/200ml



P1030001
Syringe Barrels with Pistons
Size: 3/5/10/20/30/50/60/100ml

INFUSION PRODUCT




P1008001-I
IV Cannula-injection valve
Size:14~26G




P1008001-W
IV Cannula-wings
Size:14~26G




P1008001-S
IV Cannula-small wings
Size:14~26G




P1008000-P
IV Cannula-pen like
Size:14~26G



P1007000
Spinal Needle
Size:16~27G



P1008000
Epidural Needle
Size:15-18G



P1008000
Scalp Vein Set
Size:18~27G



P1017001
Blood Lancet
Size:32*0.6*3MM



P1010001
Blood Collection Tube
SIZE : 3/4/5/6/7/10ml



P1016001
Central Venous Catheter Kit



P1016002
Haemodialysis Catheter Kit



P1018001
Combined Spinal and Epidural kit
Size: 16/18G*80mm



P1001001-LL
Disposable Infusion Set
Flashball, luer slip



P1001001-LC
Disposable Infusion Set
Y-site, luer lock



P1004001
Disposable Infusion Set
Flow regulator, luer lock



P1001004-LL
Disposable Infusion Set
Needle free connector, LC



P1002001-V
Disposable Transfusion Set
98ml/1700mm, with needle




P1006001
Infusion Set with Burette
100ml/120ml/150ml




P1014001
Dialysis Tube Set
20ml/30ml/50ml



P1007001
Extension tube/Components
all kinds components



P1012001
Disposable Infusion Bag
50/100/150/250/500/1000/2000/3000/5000ml



P1016001
AV Fistula Needles
15~17G




P1012002
Blood Bag
Single/Double/Triple/Durable etc



P1015001
Blood Dialyzer
High/Low Flux, 1.2/1.4/1.6/1.8m²


INFUSION PRODUCT




P1007001-F




P1007001-N




P1007001-T



P1007001-N



P1007002-T



P1007002-F



P1007003-A



P1007003-B



P1007003-C



P1013002



P1009009-A



P1009009-B



P1009003



P1009004



P100900-C1



P100900-C2



P1009002



P1009013



P100900-Y1



P100900-Y2



P1009010



P1009011



FOR ORDERS
Email: supplies@kumis.co.za
www.kumis.co.za