


UROLOGY





HS-C01 2-way Latex Foley Catheter

Material: Natural latex with silicone coated

Type: Plastic valve, Rubber valve

Specifications:

Fr6-Fr10 (Pediatric) Length: 270mm Balloon Capacity: 3cc/5cc

Fr12-Fr24(Standard) Length: 400mm Balloon Capacity: 5cc/10cc/30cc

Fr26-Fr30(Standard) Length: 400mm Balloon Capacity: 10cc/30cc

Color: Coded for visualization of size

Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn

Carton Size: 52x42x47.5cm



HS-C02 3-way Latex Foley Cathete

Material: Natural latex with silicone coated

Type: Plastic valve, Rubber valve

Specifications:

Length: 400mm

Fr16-Fr24(Standard) Balloon Capacity: 30cc

Fr26-Fr28(Standard) Balloon Capacity: 10cc/30cc

Color: Coded for visualization of size

Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn

Carton Size: 52x42x47.5cm



HS-C03 2-Way Female Foley Catheter

Material: Natural latex with silicone coated

Type: Plastic valve, Rubber valve

Specifications:

Length: 240mm


Fr12-Fr22 Balloon Capacity: 5cc/10cc/30cc

Fr24-Fr26 Balloon Capacity: 10cc/30cc

Color: Coded for visualization of size

Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn

Carton Size: 52x42x47.5cm



HS-C05 2-Way Tiemann Foley Catheter

Material: Natural latex with silicone coated

Type: Plastic valve, Rubber valve

Specifications:

Length: 400mm


Fr12-Fr22 Balloon Capacity: 5cc/10cc/30cc

Fr24 Balloon Capacity: 10cc/30cc

Color: Coded for visualization of size


Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn

Carton Size: 52x42x47.5cm



HS-C06 3-Way Double Balloons Foley Catheter

Material: Natural latex with silicone coated
 Type: Plastic valve, Rubber valve
 Specifications:
 Length: 400mm
 Size: Fr16-Fr24
 Balloon Capacity: 30cc(front)/50cc(rear)
 Color: Coded for visualization of size
 Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x42x47.5cm



HS-C08 4-way Double Balloons Foley Catheter

Material: Natural latex with silicone coated
 Specifications:
 Size: Fr16-Fr24 Length: 400mm
 Balloon Capacity: 30cc(front)/50cc(rear)
 Color: Coded for visualization of size
 Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x42x47.5cm




HS-C60 2-way Tiemann silicone foley catheter

Material: high quality 100% silicone, medical grade
 Type: Hard valve, With X-ray or without is available
 Specifications:
 Size: Fr14-Fr26
 Color: Coded for visualization of size.
 Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x43x39.5cm



HS-C09 Urethral Catheter(Red Latex)

Material: Natural latex with silicone coated
 Specifications: Fr12-Fr40
 Length: 400mm
 Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x43x47.5cm




HS-C07 2-way All Silicone Foley Catheter

Material: high quality 100% silicone, medical grade
 Type: Hard valve, With X-ray or without is available
 Specifications:
 Color: Coded for visualization of size.
 Fr6-Fr10(Pediatric) Length: 310mm Balloon Capacity: 3cc/5cc
 Fr12-Fr18(Standard) Length: 400mm Balloon Capacity: 5cc/10cc
 Fr20-Fr24(Standard) Length: 400mm Balloon Capacity: 5cc/10cc/30cc
 Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x43x39.5cm



HS-C10 Malecot Catheter

Material: Natural latex with silicone coated
 Specifications: Fr6-Fr24
 Length: 400mm
 Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x42x47.5cm



HS-C61 3-Way All Silicone Foley Catheter

Material: high quality 100% silicone, medical grade
 Type: Hard valve, With X-ray or without is available
 Specifications:
 Size: Fr18-Fr24
 Color: Coded for visualization of size.
 Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x43x39.5cm



HS-C11 Economic Urinary Drainage Bag Without Outlet

Capacity: 2000ml
 Without bottom outlet, with 90cm inlet tube,
 With non-return valve.
 Sterilization: EO
 Packing: 1pc/PE bag, 10pcs/bag, 250pcs/ctn
 Carton Size: 40x36x37cm



HS-C12 Economic Urinary Drainage Bag With Push Pull Valve

Capacity: 2000ml
 With bottom outlet, with Push Pull Valve,
 With 90cm inlet tube, with non-return valve.
 Sterilization: EO
 Packing: 1pc/PE bag, 250pcs/ctn
 Carton Size: 40x36x37cm



HS-C13 Economic Urinary Drainage Bag with T valve

Capacity: 2000ml
 With bottom outlet, with T valve,
 With 90cm inlet tube, with non-return valve.
 Sterilization: EO
 Packing: 1pc/PE bag, 250pcs/ctn
 Carton Size: 40x36x37cm



HS-C62 Economic Urinary Drainage Bag Push Pull Screw Valve

Capacity: 2000ml
 With bottom outlet, with screw valve,
 With 90cm inlet tube, with non-return valve.
 Sterilization: EO
 Packing: 1pc/PE bag, 250pcs/ctn
 Carton Size: 40x36x37cm



HS-C15 Luxury Urinary Drainage Bag

Capacity: 2000ml
 Type: square bottom, latex free, with pull push valve, with white sample port, with anti-reflux valve, with plastic tie
 Inlet tube: O.D 10 mm, 100cm length
 Sterilization: EO
 Packing: 1 pc/PE or Blister pack, 100pcs/ctn
 Carton Size: 54x36x32cm



HS-C17 Luxury Urinary Drainage Bag

Capacity: 2000ml
 Type: square bottom, latex free
 With T valve, with white sample port, With anti-reflux valve, without plastic tie
 Inlet tube: O.D 10 mm, 100cm length
 Sterilization: EO
 Packing: 1pc/Blister pack, 100pcs/ctn
 Carton Size: 54x36x32cm



HS-C18 Luxury Urinary Drainage Bag

Capacity: 2000ml
 Type: Square bottom, latex free
 With T valve,
 With white sample port,
 With anti-reflux valve, with plastic tie
 Inlet tube: O.D 10 mm, 100cm length
 Sterilization: EO
 Packing: 1pc/Blister pack, 100pcs/ctn
 Carton Size: 54x36x32cm



HS-C20 Luxury Urinary Drainage Bag

Capacity: 2000ml
 Type: square bottom, latex free
 With Single reversal valve, with needle sample port
 With anti-reflux valve, with double plastic hanger and single rope, with blue bed sheet clip.
 Inlet tube: O.D 10 mm, 100cm length
 Sterilization: EO
 Packing: 1pc/Blister pack, 40pcs/ctn
 Carton Size: 43x32x38cm



HS-C22 Urinary Drainage Bag

Capacity: 2000ml
 Type: Pyriform, latex free
 With T valve, with needle sample port, with anti-reflux valve,
 With double plastic hanger and rope,
 With white bed sheet clamp.
 Inlet tube: O.D 10 mm, 110cm length
 Sterilization: EO
 Packing: 1pc/Blister pack, 40pcs/ctn
 Carton Size: 43x30x38cm



HS-C23 Urinary Drainage Bag

Capacity: 2000ml
 Type: Pyriform, latex free
 With T valve, with needle sample port, with anti-reflux valve,
 with blue plastic hanger and rope,
 with white bed sheet clamp.
 Inlet tube: O.D 10 mm, 100cm length
 Sterilization: EO
 Packing: 1pc/Blister pack, 20pcs/ctn
 Carton Size: 31x26x22cm



HS-C25 Urinary Drainage Bag

Capacity: 2000ml

Type: Pyriform, latex free

With T valve, with needle sample port,

With anti-reflux valve,

With white double plastic hanger and rope,

With white bed sheet clamp.

Inlet tube: O.D 10mm, 110cm length

Sterilization: EO

Packing: 1pc/Blister pack, 40pcs/ctn

Carton Size: 43x30x38cm



HS-C26 Luxury Urinary Drainage Bag

Capacity: 4000ml

With T valve, Latex free,

With anti-reflux drip chamber(one part), with needle sample

port and tubing clip is optional, with anti-reflux valve,

With reinforced double plastic hanger and rope

With blue bed sheet clamp

Inlet tube: O.D 11mm, 110cm length

Sterilization: EO

Packing: 1pc/Blister pack, 40pcs/ctn

Carton Size: 43x30x38cm



HS-C28 Luxury Urinary Drainage Bag

Capacity: 2000ml

With T valve, latex free,

With anti-reflux drip chamber(three part)

With needle sample port and tubing clip is optional

With anti-reflux valve,

With reinforced double plastic hanger and rope


With blue bed sheet clamp

Inlet tube: O.D 11 mm, 110cm length

Sterilization: EO

Packing: 1pc/Blister pack, 40pcs/ctn

Carton Size: 43x30x38cm



HS-C29 Luxury Urinary Drainage Bag

Capacity: 2000ml

With T valve, latex free,

With anti-reflux drip chamber(three part)

With needle sample port and tubing clip is optional

With anti-reflux valve,

With reinforced double plastic hanger and rope


With blue bed sheet clamp

Inlet tube: O.D 11 mm, 110cm length

Sterilization: EO

Packing: 1pc/Blister pack, 40pcs/tn

Carton Size: 43x30x38cm



HS-C33 Urinary Meter Drainage Bag(Double Chambers)

Capacity: 2600ml + 400ml
 With T valve on urine meter
 With twist valve on urine meter
 With needle sample port, with anti-reflux valve,
 With reinforced double plastic hanger and rope,
 With white bed sheet clamp.
 Inlet tube: O.D 11 mm, 130cm length
 Sterilization: EO
 Packing: 1pc/Blister pack, 10pcs/ctn
 Carton Size: 48x41x25cm



HS-C31 Urinary Meter Drainage Bag

Capacity:
 Type A 2600ml+500ml, Type B 2000ml+200ml
 With T valve on urine meter, with twist valve on urine meter
 With needle sample port, with anti-reflux valve,
 With reinforced double plastic hanger and rope,
 With white bed sheet clamp.
 Inlet tube: O.D 10 mm, 120cm length
 Sterilization: EO
 Packing: 1pc/Blister pack, 20pcs/ctn
 Type A Carton Size:: 68x38x28cm
 Type B Carton Size::48x39x31cm



HS-C36 Leg Bag

Capacity: 500ml and 1000ml
 Type: ellipse shape
 With green twist valve and green tip cover,
 With anti-reflux valve, with double high elastic belts
 Sterilization: EO
 Packing: 1pc/PE pack, 200pcs/ctn
 Carton Size: 55x40x28cmcm



HS-C35 Leg Bag

Capacity: 600ml
 Type: Square
 With blue pull-push valve and blue tip cover
 With anti-reflux valve, with double high elastic belts
 Sterilization: EO
 Packing: 1pc/PE pack, 48pcs/ctn
 Carton Size: 30x26x15cm



HS-C37 Disposable Pediatric Urine Collector

Size: 100ml, 200ml
 Component: Pectic paper, PVC bag, sponge
 Sterilization: EO
 Packing: 1pc/PE bag, 100pcs/box, 5000pcs/ctn
 Carton Size: 55x54x43cm



FOR ORDERS
Email: supplies@kumis.co.za
www.kumis.co.za