


WOUND CARE






Wound Plaster / Adhesive Bandage

Available in different material
Available in different size and shape
H-shape/Butterfly-shape

Elastic/PU (Polyurethane)/Elastic waterproof/Perforated PE
70mm x 18mm, 72mm x 19mm 72mm x 25mm

*Other specifications or package can be produced as per requirements



Transparent PU Dressing

100 % waterproof ,thin & transparent
Adopting advanced PU material,breathable and water-resistant,
Steam permeable 10000-20000mm2/24h.

PU+Non woven Pad
Acrylic acid glue of 35-37g/M2
10dmm thickness

*Other specifications or package can be produced as per requirements



Island Wound Dressing

Sterile/Non-sterile
Paper/paper-film pouch or blister pack
Different shape and size

10cm x 8cm 10cm x 15cm 10cm x 20cm
10cm x 30cm

*Other specification or package can be produced as per requirements



Foam Dressing

Made of medical polyurethane foam with CMC , including non-woven foam dressing, foam dressing with polyurethane film, self-adherent foam dressing and silver ion foam dressing

Product thickness:2-6 mm,GSM:735g/m2/5mm(thickness)
The moisture permeability of polyurethane film which laminated on the foam dressing: 1.84g/10cm2/24h/37℃
Absorbet capability 6.94g/10cm2/24h/37℃
Total absorbent capacitance 8.78g/10cm2/24h/37℃

Hygroscopic ratio≥8



Adhesive Tape(PE/Nonwoven/Silk)

Composed of PE film/nonwoven/silk and medical pressure sensitive adhesive For the fixation of all kinds of dressings, syringe needles, catheters, etc

Transparent color
1.25cm/2.5cm / 5cm / 7.5cm / 10cm x5y
1.25cm/2.5cm / 5cm / 7.5cm/10cmx10y

Simple packing, plastic packing, metal packing

*Other specifications or package can be produced as per requirements



Perforated Zinc Oxide Adhesive Tape

Perforated zinc oxide plaster, using special processing to increase the viscosity and permeability of products

Intend usage: Protect fingers, hands, arms, knees, for the fixation of all kinds of dressings, syringe needles, catheters, etc

Base material: Cotton
Available color: White and skin

5cm / 10cm / 12cm / 18cm x5m

*Other specifications or package can be produced as per requirements



Cohesive Elastic Bandage

Easy Tear Nonwoven Fabric

Available in different color

Available in different size

1" x 5yds

2" x 5yds

3" x 5yds

4" x 5yds

6" x 5yds

*Other specifications or package can be produced as per requirements



Zinc Oxide Plaster

Premium fabric

Zinc oxide adhesive, low irritation risk

Easy hand tear, quicker use

Soft and breathable

Easy unwinds to core



Non-woven Fixing Tape

Composed of silk and medical pressure sensitive adhesive

Spunlaced non-woven tape is mainly has three different material:
common non-woven, water-repellency and waterproof

5cmx10m, 7.5cmx10m, 10cmx10m,
10cmx10m, 20cmx10m, 30cmx10m

1 Roll/Box, 150,125,80,50,25 boxes/carton

*Other specifications or package can be produced as per requirements



FOR ORDERS
Email: supplies@kumis.co.za
www.kumis.co.za