

SKIN – AESTHETIC DEVICES



Hyaluronic Acid Dermal Filler

MONALISA

Hyaluronic Acid Dermal Filler with Lidocaine

MONALISA *Lidocaine*



- Made from high molecular weight hyaluronic acid by microbial engineering
- A sterile and transparent gel of stabilized hyaluronic acid of non-animal origin
- Supplied in the pre-filled syringe with thin-wall needles

*Long-term
Stability*



Hyaluronic Acid Dermal Filler

MONALISA

MONALISA provides superior volume effect and maintains high volume for long-lasting durability.



Product	SOFT	MILD	HARD
REF	HAF10S27T	HAF10M27T	HAF10H27T
Appearance	No impurities, transparent and colorless gel		
Composition	Cross-linked hyaluronic acid		
Concentration	20 mg/mL		
Particle Size	200 µm	400 µm	600 µm
Syringe Volume	1.0 mL		
Recommended Indication	Superficial dermis	Superficial dermis / Middle layer of subcutis	Middle to deep layer of subcutis
Needle Size*	27G TW (2ea)	27G TW (2ea)	27G TW (2ea)
Storage	2~25°C		

* Needle size is subject to change without notice.

Hyaluronic Acid Dermal Filler with Lidocaine

MONALISA *Lidocaine*

MONALISA Lidocaine filler has been developed by adding local anesthetic lidocaine to lessen pain during an operation.



Product	SOFT	MILD	HARD	ULTRA
REF	MPF10S	MPF10M	MPF10H	MPF10U
Appearance	No impurities, transparent and colorless gel			
Composition	Cross-linked hyaluronic acid			
Concentration	24 mg/mL			
Lidocaine	0.3 %			
Particle Size	200 µm	400 µm	600 µm	900 µm
Syringe Volume	1.0 mL			
Recommended Indication	Superficial dermis	Superficial dermis / Middle layer of subcutis	Middle to deep layer of subcutis	Deep to very deep layer of subcutis
Needle Size*	30G TW (2ea)	27G TW (2ea)	25G TW , 27G TW	25 G TW , 27G TW
Storage	2~25°C			

* Needle size is subject to change without notice.



FOR ORDERS
Email: supplies@kumis.co.za
www.kumis.co.za